

TEFF1 Polyclonal Antibody

Cat No: HR1AP12219

For research use only

Overview

Product Name	TEFF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	TMEFF1 C9orf2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:May inhibit NODAL and BMP signaling during neural patterning (By similarity). May be a tumor suppressor in brain cancers.,similarity:Belongs to the tomoregulin family.,similarity:Belongs to the UPF0439 family.,similarity:Contains 1 EGF-like domain.,similarity:Contains 2 Kazal-like domains.,subunit:May interact with ST14.,tissue specificity:Expressed in brain.,tissue specificity:Expressed predominantly in brain, and at lower levels in heart, placenta and skeletal muscle. Down-regulated in brain tumors as compared to control brain tissues.,</p>