

ECE1 Polyclonal Antibody

Cat No: HR1AP12216

For research use only

Overview

Product Name	ECE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	ECE1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>endothelin converting enzyme 1(ECE1) Homo sapiens The protein encoded by this gene is involved in proteolytic processing of endothelin precursors to biologically active peptides. Mutations in this gene are associated with Hirschsprung disease, cardiac defects and autonomic dysfunction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.[provided by RefSeq, Sep 2009],</p>