

## SPIT2 Polyclonal Antibody

Cat No: HR1AP12209

For research use only

### Overview

Product Name	SPIT2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	SPINT2 HAI2 KOP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serine peptidase inhibitor, Kunitz type 2(SPINT2) Homo sapiens This gene encodes a transmembrane protein with two extracellular Kunitz domains that inhibits a variety of serine proteases. The protein inhibits HGF activator which prevents the formation of active hepatocyte growth factor. This gene is a putative tumor suppressor, and mutations in this gene result in congenital sodium diarrhea. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2009],