

FCRL1 Polyclonal Antibody

Cat No: HR1AP12204

For research use only

Overview

Product Name	FCRL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	FCRL1 FCRH1 IFGP1 IRTA5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Fc receptor like 1(FCRL1) Homo sapiens This gene encodes a member of the immunoglobulin receptor superfamily and is one of several Fc receptor-like glycoproteins clustered on the long arm of chromosome 1. The encoded protein contains three extracellular C2-like immunoglobulin domains, a transmembrane domain and a cytoplasmic domain with two immunoreceptor-tyrosine activation motifs. This protein may play a role in the regulation of cancer cell growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009].</p>