

UTER Polyclonal Antibody

Cat No: HR1AP12200

For research use only

Overview

Product Name	UTER Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	10kDa
GeneID?Human?	SCGB1A1 CC10 CCSP UGB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	secretoglobin family 1A member 1(SCGB1A1) Homo sapiens This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma. [provided by RefSeq, May 2010],