

VWCE Polyclonal Antibody

Cat No: HR1AP12201

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | VWCE Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 105kDa |
| GeneID?Human? | VWCE URG11 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | function:May be a regulatory element in the beta-catenin signaling pathway and a target for chemoprevention of hepatocellular carcinoma.,induction:By HBxAg.,similarity:Contains 4 EGF-like domains.,similarity:Contains 6 VWFC domains.,subcellular location:According to PubMed:16496348 is localized in the cytoplasm.,tissue specificity:Expressed in liver., |