

SEM6B Polyclonal Antibody

Cat No: HR1AP12198

For research use only

Overview

Product Name	SEM6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	97kDa
GeneID?Human?	SEMA6B SEMAN SEMAZ UNQ1907/PRO4353
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	semaphorin 6B(SEMA6B) Homo sapiens This gene encodes a member of the semaphorin family, a group of proteins characterized by the presence of a conserved semaphorin (sema) domain. Whereas some semaphorins are transmembrane proteins, others are secreted. Semaphorins play a major role in axon guidance. The protein encoded by this gene may be involved in both peripheral and central nervous system development. [provided by RefSeq, Jul 2008],