

## SIG11 Polyclonal Antibody

Cat No: HR1AP12196

For research use only

### Overview

Product Name	SIG11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	SIGLEC11 UNQ9222/PRO28718
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>sialic acid binding Ig like lectin 11(SIGLEC11) Homo sapiens This gene encodes a member of the sialic acid-binding immunoglobulin-like lectin family. These cell surface lectins are characterized by structural motifs in the immunoglobulin (Ig)-like domains and sialic acid recognition sites in the first Ig V set domain. This family member mediates anti-inflammatory and immunosuppressive signaling. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],</p>