

MMEL1 Polyclonal Antibody

Cat No: HR1AP12194

For research use only

Overview

Product Name	MMEL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	MMEL1 MELL1 MMEL2 NEP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>membrane metalloendopeptidase like 1(MMEL1) Homo sapiens The protein encoded by this gene is a member of the neutral endopeptidase (NEP) or membrane metallo-endopeptidase (MME) family. Family members play important roles in pain perception, arterial pressure regulation, phosphate metabolism and homeostasis. This protein is a type II transmembrane protein and is thought to be expressed as a secreted protein. This gene is expressed mainly in testis with weak expression in the brain, kidney, and heart. [provided by RefSeq, Jul 2008].</p>