

TMPS5 Polyclonal Antibody

Cat No: HR1AP12193

For research use only

Overview

Product Name	TMPS5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	TPRSS5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transmembrane protease, serine 5(TMPS5) Homo sapiens This gene encodes a protein that belongs to the serine protease family. Serine proteases are known to be involved in many physiological and pathological processes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014],