

# RETNB Polyclonal Antibody

Cat No: HR1AP12191

For research use only

## Overview

Product Name	RETNB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	RETNLB CCRG FIZZ2 HXCP2 RETNL2 UNQ408/PRO770
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Probable hormone.,similarity:Belongs to the resistin/FIZZ family.,subunit:Homodimer; disulfide-linked.,tissue specificity:Expressed only in the gastrointestinal tract, particularly the colon.,