

PAR12 Polyclonal Antibody

Cat No: HR1AP12190

For research use only

Overview

Product Name	PAR12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	PARP12 ZC3HDC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,similarity:Contains 2 WWE domains.,similarity:Contains 4 C3H1-type zinc fingers.,