

# KLRF1 Polyclonal Antibody

Cat No: HR1AP12176

For research use only

## Overview

Product Name	KLRF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	KLRF1 CLEC5C ML
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	killer cell lectin like receptor F1(KLRF1) Homo sapiens KLRF1, an activating homodimeric C-type lectin-like receptor (CTLR), is expressed on nearly all natural killer (NK) cells and stimulates their cytotoxicity and cytokine release (Kuttruff et al., 2009 [PubMed 18922855]).[supplied by OMIM, Oct 2009].