

NGN2 Polyclonal Antibody

Cat No: HR1AP12175

For research use only

Overview

Product Name	NGN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	NEUROG2 ATOH4 BHLHA8 NGN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neurogenin 2(NEUROG2) Homo sapiens This gene encodes a neural-specific basic helix-loop-helix (bHLH) transcription factor that can specify a neuronal fate on ectodermal cells and is expressed in neural progenitor cells within the developing central and peripheral nervous systems. The protein product of this gene also plays a role in the differentiation and survival of midbrain dopaminergic neurons. [provided by RefSeq, Apr 2012],</p>