

RSPO2 Polyclonal Antibody

Cat No: HR1AP12174

For research use only

Overview

Product Name	RSPO2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	RSPO2 UNQ9384/PRO34209
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	R-spondin 2(RSPO2) Homo sapiens This gene encodes a member of the R-spondin family of proteins. These proteins are secreted ligands of leucine-rich repeat containing G protein-coupled receptors that enhance Wnt signaling through the inhibition of ubiquitin E3 ligases. A chromosomal translocation including this locus that results in the formation of a gene fusion has been identified in multiple human cancers. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],