

CARF Polyclonal Antibody

Cat No: HR1AP12173

For research use only

Overview

Product Name	CARF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	CDKN2AIP CARF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	CDKN2A interacting protein(CDKN2AIP) Homo sapiens The protein encoded by this gene regulates the DNA damage response through several different signaling pathways. One such pathway is the p53-HDM2-p21(WAF1) pathway, which is critical to the DNA damage response. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2015],