

MXRA5 Polyclonal Antibody

Cat No: HR1AP12160

For research use only

Overview

Product Name	MXRA5 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	311kDa
GeneID?Human?	MXRA5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	matrix remodeling associated 5(MXRA5) Homo sapiens This gene encodes one of the matrix-remodelling associated proteins. This protein contains 7 leucine-rich repeats and 12 immunoglobulin-like C2-type domains related to perlecan. This gene has a pseudogene on chromosome Y. [provided by RefSeq, Mar 2010],