

## CREG1 Polyclonal Antibody

Cat No: HR1AP12158

For research use only

### Overview

Product Name	CREG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	CREG1 CREG UNQ727/PRO1409
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cellular repressor of E1A stimulated genes 1(CREG1) Homo sapiens The adenovirus E1A protein both activates and represses gene expression to promote cellular proliferation and inhibit differentiation. The protein encoded by this gene antagonizes transcriptional activation and cellular transformation by E1A. This protein shares limited sequence similarity with E1A and binds both the general transcription factor TBP and the tumor suppressor pRb in vitro. This gene may contribute to the transcriptional control of cell growth and differentiation. [provided by RefSeq, Jul 2008],