

## SEGN Polyclonal Antibody

Cat No: HR1AP12156

For research use only

### Overview

Product Name	SEGN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	SCGN SECRET
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	secretagogin, EF-hand calcium binding protein(SCGN) Homo sapiens The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation. [provided by RefSeq, Jul 2008],