

# THIOM Polyclonal Antibody

Cat No: HR1AP12153

For research use only

## Overview

Product Name	THIOM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	TXN2 TRX2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	thioredoxin 2(TXN2) Homo sapiens This nuclear gene encodes a mitochondrial member of the thioredoxin family, a group of small multifunctional redox-active proteins. The encoded protein may play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis. [provided by RefSeq, Jul 2008],