

GBX2 Polyclonal Antibody

Cat No: HR1AP12140

For research use only

Overview

Product Name	GBX2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	GBX2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May act as a transcription factor for cell pluripotency and differentiation in the embryo.,similarity:Contains 1 homeobox DNA-binding domain.,