

## DPPA4 Polyclonal Antibody

Cat No: HR1AP12139

For research use only

### Overview

Product Name	DPPA4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	DPPA4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in maintaining cell pluripotentiality.,