

GRASP Polyclonal Antibody

Cat No: HR1AP12138

For research use only

Overview

Product Name	GRASP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	GRASP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>general receptor for phosphoinositides 1 associated scaffold protein(GRASP) Homo sapiens This gene encodes a protein that functions as a molecular scaffold, linking receptors, including group 1 metabotropic glutamate receptors, to neuronal proteins. The encoded protein contains conserved domains, including a leucine zipper sequence, PDZ domain and a C-terminal PDZ-binding motif. Alternately spliced transcript variants have been observed for this gene.[provided by RefSeq, Dec 2012],</p>