

## DYN3 Polyclonal Antibody

Cat No: HR1AP12137

For research use only

### Overview

Product Name	DYN3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	DNM3 KIAA0820
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	dynamain 3(DNM3) Homo sapiens This gene encodes a member of a family of guanosine triphosphate (GTP)-binding proteins that associate with microtubules and are involved in vesicular transport. The encoded protein functions in the development of megakaryocytes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013],