

SIG10 Polyclonal Antibody

Cat No: HR1AP12136

For research use only

Overview

Product Name	SIG10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	SIGLEC10 SLG2 UNQ477/PRO940
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>sialic acid binding Ig like lectin 10(SIGLEC10) Homo sapiens SIGLECs are members of the immunoglobulin superfamily that are expressed on the cell surface. Most SIGLECs have 1 or more cytoplasmic immune receptor tyrosine-based inhibitory motifs, or ITIMs. SIGLECs are typically expressed on cells of the innate immune system, with the exception of the B-cell expressed SIGLEC6 (MIM 604405).[supplied by OMIM, Jul 2002],</p>