

BNI3L Polyclonal Antibody

Cat No: HR1AP12134

For research use only

Overview

Product Name	BNI3L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	BNIP3L BNIP3A BNIP3H NIX
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BCL2 interacting protein 3 like(BNIP3L) Homo sapiens This gene encodes a protein that belongs to the pro-apoptotic subfamily within the Bcl-2 family of proteins. The encoded protein binds to Bcl-2 and possesses the BH3 domain. The protein directly targets mitochondria and causes apoptotic changes, including loss of membrane potential and the release of cytochrome c. [provided by RefSeq, Feb 2015],