

# ISG15 Polyclonal Antibody

Cat No: HR1AP12129

For research use only

## Overview

|                       |  |
|-----------------------|--|
| Product Name          | ISG15 Polyclonal Antibody  |
| Source                | Rabbit   |
| Applications          | WB,ELISA   |
| Species Reactivity    | Human,Mouse  |
| Recommended Dilutions |  |
| Immunogen             |  |
| Species               | Rabbit   |
| Storage               | -20°C/1 year   |
| Isotype               |  |
| Clonality             |  |
| Concentration         | 1 mg/ml  |
| Observed band         | 18kDa  |
| GeneID?Human?         | ISG15 G1P2 UCRP  |
| Human Swiss-Prot No.  |  |
| Cellular localization |  |
| Alternative Names     |  |
| Background            | <p>ISG15 ubiquitin-like modifier(ISG15) Homo sapiens The protein encoded by this gene is a ubiquitin-like protein that is conjugated to intracellular target proteins upon activation by interferon-alpha and interferon-beta. Several functions have been ascribed to the encoded protein, including chemotactic activity towards neutrophils, direction of ligated target proteins to intermediate filaments, cell-to-cell signaling, and antiviral activity during viral infections. While conjugates of this protein have been found to be noncovalently attached to intermediate filaments, this protein is sometimes secreted. [provided by RefSeq, Dec 2012].</p> |