

FACE2 Polyclonal Antibody

Cat No: HR1AP12125

For research use only

Overview

Product Name	FACE2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	RCE1 FACE2 RCE1A RCE1B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Ras converting CAAX endopeptidase 1(RCE1) Homo sapiens This gene encodes an integral membrane protein which is classified as a member of the metalloproteinase family. This enzyme is thought to function in the maintenance and processing of CAAX-type prenylated proteins. [provided by RefSeq, Jul 2008],