

LRRN1 Polyclonal Antibody

Cat No: HR1AP12123

For research use only

Overview

Product Name	LRRN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	78kDa
GeneID?Human?	LRRN1 KIAA1497 Nbla10449 UNQ693/PRO1338
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 fibronectin type-III domain.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 11 LRR (leucine-rich) repeats.,