

## DPPA5 Polyclonal Antibody

Cat No: HR1AP12119

For research use only

### Overview

Product Name	DPPA5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	DPPA5 ESG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	developmental pluripotency associated 5(DPPA5) Homo sapiens This gene encodes a protein that may function in the control of cell pluripotency and early embryogenesis. Expression of this gene is a specific marker for pluripotent stem cells. Pseudogenes of this gene are located on the short arm of chromosome 10 and the long arm of chromosomes 14 and 19. [provided by RefSeq, Dec 2010],