

DNER Polyclonal Antibody

Cat No: HR1AP12116

For research use only

Overview

Product Name	DNER Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	DNER BET UNQ262/PRO299
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Activator of the NOTCH1 pathway. May mediate neuron-glia interaction during astrocytogenesis.,similarity:Contains 1 follistatin-like domain.,similarity:Contains 10 EGF-like domains.,subcellular location:Present on the membrane of dendrites and cell bodies but excluded from axonal membrane. Also found in early endosomes in the somatodendritic region.,subunit:Interacts with AP1G1. Interacts with NOTCH1.,tissue specificity:Expressed in brain, spinal cord and adrenal gland.,</p>