

CBPA2 Polyclonal Antibody

Cat No: HR1AP12114

For research use only

Overview

Product Name	CBPA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	CPA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carboxypeptidase A2(CPA2) Homo sapiens Three different forms of human pancreatic procarboxypeptidase A have been isolated. The encoded protein represents the A2 form, which is a monomeric protein with different biochemical properties from the A1 and A3 forms. The A2 form of pancreatic procarboxypeptidase acts on aromatic C-terminal residues and is a secreted protein. [provided by RefSeq, Dec 2008],</p>