

PAPP2 Polyclonal Antibody

Cat No: HR1AP12111

For research use only

Overview

Product Name	PAPP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	197kDa
GeneID?Human?	PAPPA2 PLAC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	pappalysin 2(PAPPA2) Homo sapiens This gene encodes a member of the pappalysin family of metzincin metalloproteinases. The encoded protein cleaves insulin-like growth factor-binding protein 5 and is thought to be a local regulator of insulin-like growth factor (IGF) bioavailability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2010],