

NEST Polyclonal Antibody

Cat No: HR1AP12107

For research use only

Overview

Product Name	NEST Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	178kDa
GeneID?Human?	NES Nbla00170
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nestin(NES) Homo sapiens This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells. [provided by RefSeq, Sep 2011],