

FOXC2 Polyclonal Antibody

Cat No: HR1AP12108

For research use only

Overview

Product Name	FOXC2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	FOXC2 FKHL14 MFH1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	forkhead box C2(FOXC2) Homo sapiens This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the development of mesenchymal tissues. [provided by RefSeq, Jul 2008],