

ANTR2 Polyclonal Antibody

Cat No: HR1AP12100

For research use only

Overview

Product Name	ANTR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	ANTXR2 CMG2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	anthrax toxin receptor 2(ANTXR2) Homo sapiens This gene encodes a receptor for anthrax toxin. The protein binds to collagen IV and laminin, suggesting that it may be involved in extracellular matrix adhesion. Mutations in this gene cause juvenile hyaline fibromatosis and infantile systemic hyalinosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009],