

## VASN Polyclonal Antibody

Cat No: HR1AP12099

For research use only

### Overview

Product Name	VASN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	VASN SLITL2 UNQ314/PRO357/PRO1282
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:May act as an inhibitor of TGF-beta signaling.,PTM:N-glycosylated.,similarity:Contains 1 EGF-like domain.,similarity:Contains 1 fibronectin type-III domain.,similarity:Contains 11 LRR (leucine-rich) repeats.,subunit:Interacts with TGFB1, TGFB2 and TGFB3.,tissue specificity:Expressed at highest levels in aorta, at intermediate levels in kidney and placenta and at lowest levels in brain, heart, liver, lung and skeletal muscle. Within the aorta, the strongest expression is found in the tunica media of the proximal ascending aorta, the descending thoracic aorta, the abdominal aorta and the coronary arteries. Within the kidney, expression is found in the interstitial cells.,</p>