

HTRA1 Polyclonal Antibody

Cat No: HR1AP12098

For research use only

Overview

Product Name	HTRA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	HTRA1 HTRA PRSS11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	HtrA serine peptidase 1(HTRA1) Homo sapiens This gene encodes a member of the trypsin family of serine proteases. This protein is a secreted enzyme that is proposed to regulate the availability of insulin-like growth factors (IGFs) by cleaving IGF-binding proteins. It has also been suggested to be a regulator of cell growth. Variations in the promoter region of this gene are the cause of susceptibility to age-related macular degeneration type 7. [provided by RefSeq, Jul 2008],