

LIGO2 Polyclonal Antibody

Cat No: HR1AP12097

For research use only

Overview

Product Name	LIGO2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	LINGO2 LERN3 LRRN6C UNQ9234/PRO31993
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 13 LRR (leucine-rich) repeats.,