

VWF Polyclonal Antibody

Cat No: HR1AP12095

For research use only

Overview

Product Name	VWF Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	309kDa
GeneID?Human?	VWF F8VWF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>von Willebrand factor(VWF) Homo sapiens This gene encodes a glycoprotein involved in hemostasis. The encoded preproprotein is proteolytically processed following assembly into large multimeric complexes. These complexes function in the adhesion of platelets to sites of vascular injury and the transport of various proteins in the blood. Mutations in this gene result in von Willebrand disease, an inherited bleeding disorder. An unprocessed pseudogene has been found on chromosome 22. [provided by RefSeq, Oct 2015],</p>