

PINK1 Polyclonal Antibody

Cat No: HR1AP12094

For research use only

Overview

Product Name	PINK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	PINK1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTEN induced putative kinase 1(PINK1) Homo sapiens This gene encodes a serine/threonine protein kinase that localizes to mitochondria. It is thought to protect cells from stress-induced mitochondrial dysfunction. Mutations in this gene cause one form of autosomal recessive early-onset Parkinson disease. [provided by RefSeq, Jul 2008],