

PI16 Polyclonal Antibody

Cat No: HR1AP12090

For research use only

Overview

Product Name	PI16 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	PI16 CRISP9 PSPBP PSEC0164 UNQ289/PRO328
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Putative serine protease inhibitor.,miscellaneous:May serve as a marker following prostatectomy for prostate cancer.,PTM:N-glycosylated.,similarity:Belongs to the CRISP family.,subunit:Interacts with PSP94/MSMB.,tissue specificity:Expressed in prostate, testis, ovary and intestine. Concentrates in prostate cancer patient's sera.,