

## PIGR Polyclonal Antibody

Cat No: HR1AP12085

For research use only

### Overview

Product Name	PIGR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	PIGR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>polymeric immunoglobulin receptor(PIGR) Homo sapiens This gene is a member of the immunoglobulin superfamily. The encoded poly-Ig receptor binds polymeric immunoglobulin molecules at the basolateral surface of epithelial cells; the complex is then transported across the cell to be secreted at the apical surface. A significant association was found between immunoglobulin A nephropathy and several SNPs in this gene.[provided by RefSeq, Sep 2009],</p>