

C1QT1 Polyclonal Antibody

Cat No: HR1AP12080

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | C1QT1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 30kDa |
| GeneID?Human? | C1QTNF1 CTRP1 UNQ310/PRO353 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain., |