

# EXTL3 Polyclonal Antibody

Cat No: HR1AP12078

For research use only

## Overview

Product Name	EXTL3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GeneID?Human?	EXTL3 EXTL1L EXTR1 KIAA0519
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>exostosin like glycosyltransferase 3(EXTL3) Homo sapiens This gene encodes a single-pass membrane protein which functions as a glycosyltransferase. The encoded protein catalyzes the transfer of N-acetylglucosamine to glycosaminoglycan chains. This reaction is important in heparin and heparan sulfate synthesis. Alternative splicing results in the multiple transcript variants. [provided by RefSeq, Nov 2012],</p>