

NRARP Polyclonal Antibody

Cat No: HR1AP12075

For research use only

Overview

Product Name	NRARP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	NRARP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in the formation of somites.,similarity:Belongs to the NRARP family.,similarity:Contains 2 ANK repeats.,