

TXND5 Polyclonal Antibody

Cat No: HR1AP12076

For research use only

Overview

Product Name	TXND5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	TXNDC5 TLP46 UNQ364/PRO700
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>thioredoxin domain containing 5(TXNDC5) Homo sapiens This gene encodes a protein-disulfide isomerase. Its expression is induced by hypoxia and its role may be to protect hypoxic cells from apoptosis. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring upstream MUTED (muted homolog) gene. [provided by RefSeq, Dec 2010],</p>