

OTOR Polyclonal Antibody

Cat No: HR1AP12074

For research use only

Overview

Product Name	OTOR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	OTOR FDP MIAL UNQ3054/PRO9873
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>otoraplin(OTOR) Homo sapiens This gene encodes a member of the melanoma-inhibiting activity gene family. The encoded protein is secreted via the Golgi apparatus and may function in cartilage development and maintenance. A frequent polymorphism in the translation start codon of this gene can abolish translation and may be associated with forms of deafness. [provided by RefSeq, Jul 2013],</p>