

PRVA Polyclonal Antibody

Cat No: HR1AP12068

For research use only

Overview

Product Name	PRVA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	PVALB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	parvalbumin(PVALB) Homo sapiens The protein encoded by this gene is a high affinity calcium ion-binding protein that is structurally and functionally similar to calmodulin and troponin C. The encoded protein is thought to be involved in muscle relaxation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015],