

PGRP4 Polyclonal Antibody

Cat No: HR1AP12063

For research use only

Overview

Product Name	PGRP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	PGLYRP4 PGRPIB SBBI67
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>peptidoglycan recognition protein 4(PGLYRP4) Homo sapiens Summary: This gene encodes a peptidoglycan recognition protein, which belongs to the N-acetylmuramoyl-L-alanine amidase 2 family. These proteins are part of the innate immune system and recognize peptidoglycan, a ubiquitous component of bacterial cell walls. This antimicrobial protein binds to murein peptidoglycans of Gram-positive bacteria. [provided by RefSeq, Oct 2014],</p>